

SYNTHdetect



Overview

SYNTHdetect is a desktop diamond verification instrument for loose and mounted stones, designed to separate diamonds from synthetic diamonds. SYNTHdetect is manually operated and requires an operator to interpret the test results.

SYNTHdetect captures images and based on visual analysis of the images the instrument automatically displays one of the following results that can be interpreted by the operator:

i) Natural ii) Refer

The 'Refer result could contain potential synthetic diamonds, diamonds and diamond simulants. These referrals should be tested further for full verification. The instrument does not require user interpretation although it contains the software to aid in this interpretation. SYNTHdetect needs to be calibrated periodically.

Stone Testing Capabilities

Weight range:	0.001 ct - 100 ct
Colour range:	D-J failsafe across all colours although natural
	fancy colours will give a higher referral rate
Stone shape:	All
Diamond simulants:	Yes
Mounted stones:	Yes not yet ASSURE tested

Instrument Capabilities

Automated feed:	No
Automated results:	No
Automated dispense:	No
Detect or Refer (synthetic diamonds):	Refer
Detect or Refer (diamond simulants):	Refer
Multiple stones at once:	Yes
Training required:	Yes

Instrument Performance

The results demonstrate the performance of the SYNTH detect tested by an expert operator:

	Core Sample	Smalls Sample
Diamond False Positive Rate Optimal rate 0%	0.0%	0.0%
Diamond Referral Rate Optimal rate 0%	0.7%	1.9%
Diamond Accuracy Rate Optimal rate 100%	99.3%	98.1%
Speed Stones tested per hour	583	583

Instrument Specifications

31 cm (W) x 34 cm (D) x 45.8 cm (H)

- 30 kg
- Price: 18,500 USD July 2019

Contact

debeersgroupservices.com/instruments contact@debeersgroupservices.com +44 207 858 7887

ASSURE Tested March 2019 ASSURE ID 375648